

Novel microphone to facilitate communication during non-invasive ventilation (NIV)

ReddyPort™ non-invasive ventilation products

Chakravarthy B. Reddy MD, Dept. of Pulmonary and Critical Care Medicine, University of Utah

Elizabeth D. Young BS, Dept. of Communication Sciences and Disorders, University of Utah

Lara M. Brewer Ph.D., Dept. of Anesthesiology, University of Utah

Benjamin F. Schiedermaier SLP, Dept. of Communication Sciences and Disorders, University of Utah

Sarah H. Ferguson Ph.D., CCC-A, Dept. of Communication Sciences and Disorders, University of Utah

Study performed at Speech Acoustics and Perception Laboratory, Department of Communication Sciences and Disorders, University of Utah, 390 South 1530 East, Room 1201, Salt Lake City, UT 84112.

Introduction

The use of non-invasive ventilation (NIV) in hospitalized patients is expanding and it is often the first therapeutic option for patients with respiratory distress¹. A major drawback of NIV is impaired communication,^{2,3} which can result in high anxiety for the patients and the caregivers.⁴ The mask, the noise of the ventilator, and the high rate of airflow reduce intelligibility of speech. This also reduces the ability to engage in conversations around patients' goals of care. Intermittent removal of the mask can facilitate speech and communication. However, removal of the mask leads to loss of positive pressure and derecruitment of alveoli, which often causes a deterioration in respiratory status. Communication impairment assumes even more significance in patients requiring NIV for palliative support at the end of life. Innovations to improve communication in these patients are much needed.⁵

For this study, we tested a new microphone (*ReddyPort™ Microphone; ReddyPort Inc., Salt Lake City, UT*) that amplifies speech during NIV. To use the microphone, the elbow of an NIV mask is replaced with an adapter that has a one-way, self-sealing valve. This design maintains positive pressure and allows access to the oral cavity without causing any air to leak. The microphone fits into the adapter without interrupting NIV airflow or positive pressure (*figure 1*). The microphone is a single-patient-use device with a reusable controller. The microphone works by amplifying the patient's voice with an integrated speaker, while clarifying the output by using digital signal processing to remove background breathing and ventilator noise. The digital signal processing also naturalizes the patient's voice.

Funding Source

The study was funded by NIH NINR SBIR grant (FWA00024352)

Corresponding author

Chakravarthy Reddy, MD
Dept. of Pulmonary and Critical Care Medicine
University of Utah Health Sciences Center
1950 Circle of Hope
Salt Lake City, UT 84112-5550
Telephone: 801.587.4762
Email: c.reddy@hci.utah.edu

Key words

Non-invasive ventilation (NIV), communication impairment, intelligibility of speech, positive airway pressure, patient communication.



Figure 1: ReddyPort Microphone attached to an NIV mask.

Continued on the next page.

In this study, we tested the improvement in intelligibility of speech using ReddyPort Microphone in healthy volunteers on NIV. In a sound-treated booth, five healthy adult volunteers were recorded reading ten standard IEEE sentences while wearing an oronasal NIV mask (*Phillips Respironics AF541*). Study subjects were closely monitored for effects of hyperventilation. Audio recordings were obtained with the microphone off and on during three separate NIV conditions: Baseline (*no airflow*), Continuous Positive Airway Pressure (CPAP) 5 cm H₂O and Bi-level Positive Airway Pressure (BiPAP) 8/4 cm H₂O. Recorded sentences were then presented to 12 blinded listeners who transcribed what they heard. Each response was scored as a percentage of keywords correct. Linear mixed-effects models were used to analyze the effects of NIV condition and microphone setting. Between the microphone off and on conditions, sentence keyword intelligibility increased significantly from 78.6% to 85.8% with CPAP ($z = 5.16, p < 0.001$) and from 55.8% to 79.6% with BiPAP ($z = 13.41, p < 0.001$). As expected, the baseline condition with no airflow showed little change in intelligibility, from 93.5% to 89.2% ($z = -4.83, p < 0.001$) (*figure 2*). The study was approved by the Institutional Review Board of University of Utah (*IRB Nos. 56031 and 56030*). It was conducted in the Voice, Airway, and Swallowing Transitional Laboratory and the perceptual experiment was carried out by personnel from the Speech Acoustics and Perception Laboratory.

Results of intelligibility testing




Figure 2: results of intelligibility testing.

Two other devices that can improve communication while on NIV are available and have been reported.^{6,7} First is a speech enhancement device utilizing a throat sensor that converts vocal cord vibrations into an acoustic voice signal.⁶ The second is a Positive Airway Pressure communication device that was shown to improve intelligibility of speech in an out-patient setting in patients with obstructive sleep apnea using continuous positive airway pressure mask.⁷ The advantage with the microphone evaluated in this study is the opportunity to detect and adjust for in-mask ventilator noise, weak patient voices, and unnatural voice quality. Moreover, with this microphone, the improvement appears to be more pronounced with increasing pressure settings because of decreasing baseline intelligibility from the ventilator noise.

Conclusion

Using the ReddyPort Microphone integrated into the NIV mask improved sentence keyword intelligibility in healthy volunteers on NIV. Utility of this microphone in patients in an ICU setting and with higher levels of BiPAP is currently being evaluated. The microphone has the potential to provide a significant improvement towards enhancing communication in critically ill patients treated with NIV.

For more information call **801.899.3036**.

References:

- 1 Aboussouan LS, Ricarte B. Noninvasive positive pressure ventilation: Increasing use in acute care. *Cleve Clin J Med.* 2010 May;77(5):307-16.
- 2 Green E, Bernoth M. The experiences of nurses using noninvasive ventilation: An integrative review of the literature. *Aust Crit Care.* 2020 Nov;33(6):560-566.
- 3 Sørensen D, Frederiksen K, Groefte T, Lamborg K. Nurse-patient collaboration: a grounded theory study of patients with chronic obstructive pulmonary disease on non-invasive ventilation. *Int J Nurs Stud.* 2013 Jan;50(1):26-33.
- 4 Schmidt M, Boutmy-Deslandes E, Perbet S, Mongardon N, Dres M, Razazi K, Guerot E, Terzi N, Andrivet P, Alves M, Sonneviller R, Cracco C, Peigne V, Collet F, Sztrymf B, Rafat C, Reuter D, Fabre X, Labbe V, Tachon G, Minet C, Consell M, Azoulay E, Similowski T, Demoule A. Differential Perceptions of Noninvasive Ventilation in Intensive Care among Medical Caregivers, Patients, and Their Relatives: A Multicenter Prospective Study-The PARVENIR Study. *Anesthesiology.* 2016 Jun;124(6):1347-59.
- 5 Wong AI, Cheung PC, Hopp MB, Gay PC, Collop NA. Consequences and Solutions for the Impact of Communication Impairment on Noninvasive Ventilation Therapy for Acute Respiratory Failure: A Focused Review. *Crit Care Explor.* 2020 Jun 15;2(6):e0121.
- 6 Jssennaggar CE, Ten Hoorn S, Girbes AR, Tuinman PR. A new speech enhancement device for critically ill patients with communication problems: a prospective feasibility study. *Intensive Care Med.* 2017 Mar;43(3):460-462.
- 7 Wong AI, Cheung PC, Zhang J, Cotsonis G, Kutner M, Gay PC, Collop NA. Randomized Controlled Trial of a Novel Communication Device Assessed During Noninvasive Ventilation Therapy. *Chest.* 2021 Apr;159(4):1531-1539.

ReddyPort, Salt Lake City, UT, 84101, U.S.

reddyport.com

© 2021 ReddyPort and the ReddyPort logo are trademarks of ReddyPort and Company. The Philips and Respironics trademarks are property of Koninklijke Philips N.V. or their respective owners. RP5002 (0921)

REDDYPORT